



Chemical, Physical and Thermodynamic Evaluation of Arsenic Species and Their Transformation during Cigarette Smoking

C. Liu, J. Hu, C. Wright
K. McAdam

British American Tobacco
Group Research and Development,
Southampton
SO15 8TL, UK

H. Goenaga Infante, S.
Taebunpakul, J. Braybrook

LGC Ltd.
Queens Road, Teddington
Middlesex
TW11 0LY, UK

Background & Objectives

- Trace metals in smoke: important?
 - As, Cr, Cd, Ni, Pb, Se
 - Toxicity depending on: level and speciation
- Trace metals analysis: challenges
 - Matrix: over 5600 identified constituents in tobacco/smoke
 - Routine metal analysis in smoke: total level
 - Smoke collection: special care needed
- X-ray absorption spectroscopy: synchrotron radiation
- Objective: further understanding of As speciation via chemical analysis

Arsenic Toxicity: Valence & Speciation

- Inorganic:



- Arsenite (+3) LD_{50} (mg/kg) = 15 to 42
- Arsenate (+5) LD_{50} (mg/kg) = 20 to 200

- Organic: less toxic

- Monomethyl arsenic acid (MMA): +3 & +5 LD_{50} (mg/kg) = 700 to 1800
- Dimethylarsenic acid (DMA): +3 & +5 LD_{50} (mg/kg) = 1,200 to 2,600
- Arsenocholine LD_{50} (mg/kg) = 6,500
- Arsenobetaine (AsB) LD_{50} (mg/kg) \geq 10,000

(C.O. Abernathy, HECD/OST, US EPA)

Time Scale in Human Body: Arsenic



26.4 ppm in smoke (1933), 14.4 ppm in smoke (1959) (Wynder & Hoffmann: Tobacco & Tobacco Smoke), ~1.0 ppm (2008)

R.H. Holland, et al, Cancer Research, 19 (1959), pp. 1154 - 1157

SAMPLE PREPARATIONS

- Tobacco & ash preparation
 - Tobacco: cut rag from as-received 3R4F cigarette (x20)
 - Ash: ISO, machine-smoked 3R4F cigarette (x10)
- Smoke sample preparation
 - ISO, machine-smoking, metal-free substrate*
 - 'Fresh': solid-CO₂ cooled during smoking/storage/transfer
 - 'Aged': 1 day to 7 days, ambient temperature during storage/transfer
- Two dimensional HPLC with ICP-MS

*C. Liu, et al., *Spectrochim Acta B*, 2009, **64**, 1294

CORESTA CONGRESS 2008, SSPT09

C. Liu *et al.*, Shanghai, China

“X- Ray Absorption Spectroscopy for the Determination of Oxidation States of Trace Metals in Tobacco, Cigarette Smoke and Ash”

- XANES analysis of cut tobacco and ash samples suggest that As(V) was the dominant species.
- In MS smoke TPM, both As(III) and As(V) were detected. Ageing of the TPM samples promotes As(III) → As(V) transformation.

However

No species identification was possible


WATER-SOLUBLE ARSENIC IN CUT 3R4F TOBACCO: EXTRACTION

> Total As (dry weight) in homogenate of 3R4F tobacco: **$318 \pm 9 \text{ ng g}^{-1}$**

Extraction	As (ng g^{-1})	Efficiency (%)
Water	135 ± 5	42
Driselase	40 ± 1	13
SDS	25 ± 2	8
Total	200 ± 10	63

Optimisation of Extraction

- Compromise between total As in the extract and species preservation
- Comparison of 1% formic acid, water, ammonium acetate buffer, water-methanol (**optimal: water**)
- To enhance As ICP-MS transient signal intensity, optimisation of sample size to extractant ratio (**optimal: 0.2 g solid in 2 g water**)
- Optimisation of the extraction time (**optimal: 10 min by MAE**)





SEPARATION OF WATER-SOLUBLE ARSENIC SPECIES IN CUT 3R4F TOBACCO

Peak ID:


- 1: Arsenobetaine (AsB)
- 2: Arsenite
- 3: Dimethylarsenic Acid (DMA)
- 4: Monomethylarsenic acid (MMA)
- 5: Arsenate

Agreement with XANES



SMOKE CONDENSATES WATER EXTRACT

(A)




← **As(III) + As(V): 51%**
U: 41%

- 1: AsB
- 2: Arsenite: As(III)
- 3: DMA
- 4: MMA
- 5: Arsenate: As(V)
- 6: Unknown

Total As, water extract: 0.11-0.17 mg/kg

(B)





← **As(V)**

Both Peak 2 As(III) and Peak 6 are temperature sensitive.

Retention time / min


Arsenic Redox Reactions During Cigarette Smoke Formation & Ageing

“Smoke pH” and redox potential: important parameters which decide the direction of arsenic species transformation



Thermodynamically: Pourbaix Diagram of Arsenic

Smoke tar: $E = +0.17$ to $+0.24\text{V (S.C.E)} = -0.07$ to 0.00 V(S.H.E)



SUMMARY

- ◆ Sequential extraction of As from 3R4F cut tobacco achieved overall efficiency of 63%.
- ◆ Microwave assisted extraction enhanced As available from cut tobacco in water, reducing extraction time significantly (2 h to 10 min).
- ◆ Physical (XANES) and chemical (HPLC-ICP-MS) methods agree that most As in cut tobacco is As(V). **However**, the latter provides additional information on As species & distribution.
- ◆ Complementary results were also achieved on mainstream smoke condensate, pointing to a As(III) → As(V) upon ageing and temperature increase.
- ◆ Further work is ongoing to identify the remaining unknown As species in smoke (most likely, thioarsenite).

This presentation is available on bat-science.com

Thank You & Any Question

www.bat-science.com



**BRITISH AMERICAN
TOBACCO**