

EASTMAN

Improved Pressure Drop & Ventilation Instrument with Variable Flow Rates

Steven A. Wilson
Eastman Chemical Company
TSRC 2011

Objective

- Replace pressure drop and ventilation (PDV) instrument.
- Develop an instrument that can operate over a range of flow rates.
- Utilize the up-to-date low flow control systems and measurement technology.


Why the replacement?

- Need to upgrade testing equipment
- Need to measure tip pressure drops & filter ventilation at various flow rates/puff volumes
 - Various smoke testing standards
 - Increased regulatory environment
 - Tighter quality control

Comparison of smoking conditions

	Puff volume (ml)	Puff duration (sec)	Puff interval (sec)	Ventilation blockage (%)	Mean flow rate (ml/s)
ISO	35	2	60	0	17.5
Massachusetts Intense	45	2	30	50	22.5
Canadian Intense	55	2	30	100	27.5

Previous research instrument




Features of previous PDV instrument

- Design and built in 1980
- Flow rate = 17.5 ml/min
- Flow adjusted with a micro metering valve and controlled by measuring the pressure drop through a Meriam laminar flow calibration standard
- Pressure drop measured by a Validyne differential pressure cell sensor
- Ventilation air-flow measured by a Kurz air-flow meter
- Dot matrix printer



Basic pressure drop instrument diagram



Comparison of flow controllers

- Critical flow orifice
 - Non adjustable
 - Inaccurate
- Needle valve
 - Allows adjustments
 - Need to calibrate
- Flow Controllers
 - Pneumatic vs. electronic control
 - Flow rate check versus calibration
 - Environmental influences
 - Temperature
 - Relative humidity
 - Atmospheric pressure

Custom Electronic Systems, Inc. (CES)

- Corporate expertise
 - Specializing in low air flow measurement and control
 - In business for over 30 years
 - Located in Rural Hall, North Carolina

- CES has previously developed research instruments for the tobacco industry
 - Model CES 209 Flow & Totalizer
 - Model CES 228 Filter Size Controller
 - Model CES 237 Flow & Volume Indicator
 - Model CES 280 Flow Controller

Features of new PDV instrument

- Designed and built in 2010
- Flow rate of up to 100 ml/sec
- Pressure drop range of 0 – 1000 mm H₂O
- Ventilation range of 0 – 100 %
- Flow rates measured with two laminar flow elements
- Microprocessors allowing real time measurement of air flow rates
- Auto zero compensation for room temperature, relative humidity, barometric pressure, and vacuum levels
- Various circumferences possible (by changing seals)

CES 508 Instrument



Sample holder




Previous Holder




New Holder




Design of sample holder




CES 508 - front panel



CES 508 - inside sections




CES 508 - back panel




Ventilation series experiment

- Objective: To evaluate the influence of flow rate on the ventilation level.
- Conditions: Compare various flow rates.
- Samples: Cigarettes with a range of ventilations (0, 15, 30, 45, 60, & 75%).
- Does the ventilation level change with air flow?
 - If one has 45% ventilation level at 35 ml puff (ISO), what level would be observed with a 55 ml puff?

Pressure drop comparison



Ventilation comparison



Precision comparison - r & R for ranges

	Cigarette Pressure Drop Repeatability (r)	Cigarette Pressure Drop Reproducibility (R)	Cigarette Ventilation Repeatability (r)	Cigarette Ventilation Reproducibility (R)
CORESTA Method No. 4	3 mm WG	6 mm WG	-	-
CORESTA Method No. 6	1.86 – 4.33 mm WG	3.99 – 8.33 mm WG	0.86 – 1.91 %	1.88 - 2.89 %
CORESTA Method No. 41	2.3 mm WG	5.8 mm WG	-	-
CES 508 Instrument	1.4 mm WG	-	1.7 %	-

Summary

- A new pressure drop and ventilation instrument with adjustable flow rates has been designed and built.
- Various new design features were utilized to allow improved accuracy and precision.
 - Redesigned air paths ways to allows higher air flow rates
 - Microprocessors allowing real time measurement of air flow rates
 - Separated zones to reduce influence of electronics heating on air measurements
 - Compensation for room temperature, relative humidity, and barometric pressure
- Experimental results show good precision.
- A demonstration experiment shows that a filter's pressure drop and ventilation are significantly influenced by the air flow rates.

BLANK

EXTRA SLIDES

Contact Information

Custom Electronic Systems, Inc. (CES)

455 Forum Parkway

Rural Hall, North Carolina 27045

336-969-2411


<http://www.drafrange.com/>

ces@drafrange.com

Precision comparison

- CORESTA Guide No. 4
 - A User Guideline for Measurement of Pressure Drop of Cigarettes and Cigarette Filter Rods
- CORESTA Recommended Method No. 6
 - Determination of Ventilation – Definitions and Measurement Principles
- CORESTA Recommended Method No. 41
 - Determination of the Draw resistance of Cigarettes and Filter Rods

Ventilation types



References

- “The effect of cigarette design on content of phenols in mainstream tobacco smoke”, Beitr. Tabakforsch. Int., 24:187-193, 2011, by S. Dagnon, et al.
- “The Effect of Filter Ventilation on the Yield and Composition of Mainstream and Side Stream Smokes”, Beitr. Tabakforsch. Int., 10:81-90, 1980, C. L. Browne, et al.
- “Parameters Affecting the Selective Filtration of Certain Tobacco Smoke Components”, Beitr. Tabakforsch. Int., 8:145-9, 1975, by G. P. Morie, et al.
- “The Effect of Perforated Tipping Paper on the Yield of Various Smoke Components”, Beitr. Tabakforsch. Int., 7:282-7, 1974, by V. Norman.

- “The influence of water on the selective filtration of phenol in upstream versus downstream segments”, Watts & Wilson, TSRC 2008.
- “Smoke composition changes resulting from filter ventilation”, S. A. Wilson, TSRC 2001.