

Highly time-resolved two-dimensional mapping of the molecular combustion and pyrolysis product concentrations during a puff in a burning cigarette

R. Zimmermann¹, R. Hertz-Schünemann¹, S. Ehlert^{1,3}, A. Walte³, C. Liu²,
K. McAdam², S. Coburn², T. Streibel¹


*1 Joint Mass Spectrometry Centre
University of Rostock, Chair of Analytical Chemistry,
Inst. of Chemistry, Germany and Helmholtz-Zentrum
München, KMS-Cooperation Group Analysis of
Complex Molecular Systems, Germany
ralf.zimmermann@helmholtz-muenchen.de*

*2 Photonion GmbH
Schwerin, Germany
andreas.walte@photonion.de*

*3 BAT Ltd.
Southampton, UK*






Introduction: Photo ionization - SPI and REMPI



2014_TSRC41_Zimmermann.pdf

analysis
e applied:



a) V
- i
- S
- p
b) R
- io
- hi
- p

➔ Mass analysis by **time-of-flight (TOFMS)**, ion trap (ITMS) etc.

2014_TSRC41_Zimmermann.pdf

Why using (vacuum-) photo ionization for fast on-line monitoring?

- + **soft ionization** (fragmentation-free ionization of many tobacco smoke chemicals and reduced fragmentation of labile compounds)
- + **SPI - universal ionization** (incl. alkanes), IE-threshold selectivity
- + **REMPI: selective for aromatic** compounds
- + two easily **switchable** laser PI techniques: REMPI & SPI (+EI)
- + **sensitive and fast**: ppb concent. in sub-second time resolution
- + **no matrix effects**: “physical ionization”

Introduction: Photo ionization - SPI and REMPI








Photo ionization – mass spectrometry: On-line tobacco smoke analysis





Adam et al., Mitschke et al., Streibel et al.

2R4F - Mainstream




2R4F - Sidestream



480
420
360
300
240
180
120
60
0

Time [s]



68th TSRC2014(68) - Document not peer-reviewed

Photo ionization – mass spectrometry: On-line tobacco smoke analysis




Photo ionization MS smoke profiler for Industrial applications (Borgwaldt KC)

- Integrated smoke machine – SPI-TOFMS system (programmable profiles: ISO, intense, human, e-cig.) for industrial routine & research
- Quantitative puff-resolved on-line analysis of tobacco smoke toxicants (e.g., isoprene, butadiene, benzene, acetaldehyde)
- Analysis of e-cigarette and heat-not-burn product smoke composition
- Easy to use integrated system (smoke machine)
- Compact oa-TOFMS mass analyzer (R~900)
- Different VUV-lamps for SPI ion source
- EI ion source

Stand-alone photo ionization MS for flexible applications (Photonion GmbH)


- Flexible, ultra-high performance PI-TOFMS system for research and industrial process analysis
- System can be coupled to:
 - Thermal analysis,
 - **μ -probe sampling device,**
 - Smoking machine,
 - Human smoking analyzer (Gas/PM)
 - Smoke particle matter analyzer
 - Fast gas chromatography pre-separation (puff- resolved)
- High flexibility and sensitivity
- Direct Reflectron TOF mass analyzer (R~2000)
- REMPI, laser-SPI and VUV lamp SPI source
- EI ion source and switchable ionization
- Multiple coupling options.

Application of PI-MS for analysis of pyrolysis & combustion products within the cigarette



2014_TSRC41_Zimmermann.pdf


- Research **Cigarettes**: Well **defined model-system** for solid state combustion
- Tobacco: Bio mass consisting mainly of **Celluloses** (“glucose-polymer”) and **Lignin** (“phenolic-polymer”) with **alkaloid** content(nicotine)
- Very **reproducible** “automated” **machine smoking** (ISO-testing): One bell-shaped puff (35 ml/2s) per minute





P = Pyrolysis/
Distillation zone
C = Combustion zone

TSRC2014(68) - Document not peer-reviewed

Application of PI-MS for analysis of pyrolysis & combustion products within the cigarette




2014_TSRC41_Zimmermann.pdf






2014_TSRC41_Zimmermann.pdf - Document not peer-reviewed

Application of PI-MS for analysis of pyrolysis & combustion products within the cigarette




2014_TSRC41_Zimmermann.pdf




2014_TSRC41_Zimmermann.pdf
C2014(68) - Document not peer-reviewed

Application of PI-MS for analysis of pyrolysis & combustion products within the cigarette




Rich pyrolysis/combustion product spectrum: Masses up to 300 m/z detectable




SPI: Universal ionization:

→ For example: Higher molecular weight signature of fatty acids ★


mass 228 Tetradecanoic acid




mass 242 Pentadecanoic acid




mass 256 Hexadecanoic acid
(Palmitic acid)



mass 270 Heptadecanoic acid




Application of PI-MS for analysis of pyrolysis & combustion products within the cigarette




2014_TSRC41_Zimmermann.pdf

Rich pyrolysis/combustion product spectrum: Aromatics and Polycyclic Aromatic Hydrocarbons (PAH) with masses up to 500 m/z detectable




2014_TSRC41_Zimmermann.pdf - Document not peer-reviewed

Application of PI-MS for analysis of pyrolysis & combustion products within the cigarette



SPI-TOFMS


REMPI-TOFMS





2014_TSRC41_Zimmermann.pdf

TSRC2014(68) - Document not peer-reviewed

PI-MS for time- and space-resolved chemical mapping of the pyrolysis/combustion zone





Mapping: Addressing multiple position in cigarette during a puff



Pressure and temperature changes during the measurement

→ SPI Signals are normalized with respect to the nitrogen LEI signals

PI-MS for time- and space-resolved chemical mapping of the pyrolysis/combustion zone




Substance	Quantification
NO	External calibration
Benzene	External calibration
Toluene	External calibration
Acetaldehyde	Cross section: 0.20*
Butadiene	Cross section: 0.60*



Substance	Relative cross section
Acetone	Cross section: 0.31*
Isoprene	Cross section: 0.59*
Ammonia	-
Phenol	-
Nicotine	-

* Relative to benzene,


PI-MS for time- and space-resolved chemical mapping of the pyrolysis/combustion zone




5 measurements averaged value at puff end for this sampling point:
~ 50110 ppm NO and 1 ppm benzene





PI-MS for time- and space-resolved chemical mapping: Distribution maps (middle of puff)






PI-MS for time- and space-resolved chemical mapping: Distribution maps (time-resolved)





PI-MS for time- and space-resolved chemical mapping: Distribution map “movie”



- 0.4 s
c1ccccc1
(ppm)



NO
(ppm)




O₂
(%)

PI-MS for time- and space-resolved chemical mapping: NO_x formation and destruction


NO_x-Formation in bio mass (tobacco) combustion:

- Temperature not sufficient for **Thermal NO_x-formation** (*Zeldovich*-mech.)
- Rate of **Prompt NO_x-formation** (*Fenimore*-mech.) small
- Fuel NO_x-Formation** is dominating: **Heterogeneous** (from N-containing char) and **Homogeneous** (via N-volatiles as HCN etc.) fuel NO_x-Formation

PI-MS for time- and space-resolved chemical mapping: 3D-Distribution maps (middle of puff)





2014_TSRC41_Zimmermann.pdf





TSRC2014(68) - Document not peer-reviewed

PI-MS for time- and space-resolved chemical mapping: Effective formation/destruction kinetics





Effective homogeneous formation/destruction kinetics of NO and benzene in a volume increment at puff peak along center line









PI-MS for time- and space-resolved chemical mapping: Other compounds



Summary



- PI-MS (REMPI, SPI): Soft ionization method for analysis of trace gases and combustion effluents
- Puff-resolved monitoring of gaseous cigarette smoke compounds (health risk estimation tobacco research)
- PI-MS with μ -probe sampling: Analyzing the composition of the pyrolysis and combustion zones: Cigarette as “solid fuel combustion model”
- “Machine smoking” allows repetitive highly dynamic PI-MS measurements: Mapping approach \rightarrow NO_x-formation
- Spin-off company Photonion GmbH: Custom-build and standard PI-MS Systems for Industry and Research



Multi-component trace gas analysers



On-line photoionisation mass spectrometry



Dr. S. Ehlert



Dr. M. Bente




Dr. M. Saraji

Instrumental Solutions and Applications



- Flexible Systems for Research and Industry: Photo-TOF (SPI/REMPI/EI)
- On-line Analysis of Cigarette Smoke (OEM for Borgwaldt KC)
- Evolved Gas Analysis in TA (OEM for Netzsch GmbH)

Contact: walte@photonion.de



Outlook – Further microprobe-PI-MS applications: Flavor formation in coffee roasting





Sampling within coffee-bean during



Journal of MASS SPECTROMETRY

December 2013 • Volume 48 • Issue No. 12
Pages 1253–1366 • ISSN 1076–5174

wileyonlinelibrary.com/journal/jms



SPECIAL FEATURE: PERSPECTIVE

On-line process monitoring of coffee roasting by resonant laser ionisation time-of-flight mass spectrometry: bridging the gap from industrial batch roasting to flavour formation inside an individual coffee bean

By R. Hertz-Schünemann, R. Dorfner, C. Yerezian, T. Streibel and R. Zimmermann

WILEY

NEW!

EXPRESS PUBLICATION 4 WEEKS
Acceptance to online paginated publication

Tb201311p01 - Document not peer-reviewed



Dr. T. Streibel



Dr. R. Hertz-Schünemann



S. Ehlert

M. Schaffer E. Karg
T. Gröger B. Gruber
R. Hertz M. Kürsten
C. Busch M. Jennerwein
J. Hölzer F. Li
K. Mthembu L. Müller
C. Schepler A. Rheda
A. Sutherland X. Wu
S. Ly-Verdu S. Scholtes
S. Wohlfahrt C. Radischat
M. Fischer M. Antretter
C. Rüger E. Müller
T. Schwemer K. Lau
C.-G. Dragan
H. Czech



Joint Mass Spectrometry Centre

Funding:

- Uni Rostock und Helmholtz Zentrum München
- Deutsche Forschungsgemeinschaft (DFG)
- Bundesministerium für Bildung & Forschung (BMBF)
- Bayerische Forschungsförderung (BFS)
- Deutsche Gesetzl. Unfallversich. (DGUV)
- Bundeskriminalamt (BKA)
- Companies (SASOL Ltd., Netzsch GmbH, Photonion GmbH, Airsense GmbH, Borgwaldt KC etc.)
- Helmholtz-Impulse and Networking- Fonds (Virtual Helmholtz Institute - HICE)