

E-Cigarette Flavour Transfer Screening Method by GC/MS


Stuart Martin, Malcolm Saxton, [Desmond Dalton](#)

CORESTA CONGRESS – October 2018

Analytical Development Centre, British American Tobacco,
Southampton, UK


Introduction

- Importance of assessing flavour transfer
- Method overview
- Validation and results
- The Limit Test
- Summary & Conclusions
- Acknowledgements
- Questions



Importance of Assessing Flavour Transfer


- Increasing requirement from a Stewardship perspective with regards to the transfer rate of certain flavour compounds from e-liquid to e-aerosol [1].
- Product development.
- Traditionally requires method development and validation of quantitative assays.
- The main objectives of this project were:
 - To provide a cost effective and fast turnaround (< 3 weeks) screening measurement for flavour transfer from e-liquids to e-aerosols.
 - To demonstrate whether transfer rates exceed BAT Stewardship Limits for defined compounds (Limit of Concern – LOC).
- The solution: a liquid to aerosol comparison method with Limit Test.



[1] Costigan, S. & Meredith, C: An approach to ingredient screening and toxicological risk assessment of flavours in e-liquids, Regulatory Toxicology and Pharmacology, Volume 72, Issue 2, July 2015, Pages 361-369.

Method Overview


Collecting the aerosol



- Extraction solution – ethanol containing internal standard (anisole-d8 @ 1µg/mL)
- 20mL / impinger
- All extractions use same batch of extraction solution
- 80mL puff volume, 3 sec puff duration, 30 sec puff interval

Method Overview

Generating the samples and calculating transfer



Validation and Results

- Target analytes: Limonene, Neomenthol, Menthol (peppermint oil e-liquid)


Criterion	n	Description
Extraction time optimisation	3	25 puffs 20mg/puff device (peppermint oil) onto CFP, extract in 20mL ethanol, aliquots taken 5, 10, 15, 20 minutes. Analyse in SCAN and SIM modes. Determine optimum extraction time.
Sintered vs Non-sintered impingers + distribution	5	Check for differences in trapping efficiency (e.g. more or less carry-over to 2nd impinger). Target $\leq 10\%$ carry-over to back impinger.
Injection precision	5	Repeat injection of the same diluted e-liquid sample. Evaluate precision ($\leq 1\%$).
Instrument carry-over	N/A	Instrument blanks analysed between samples, compound carry-over $\leq 1\%$ assessed by peak area.
Method repeatability	5	RSD $\leq 15\%$ for all analytes



Validation and Results

- Extraction time optimisation

- 15 minutes @ 180rpm selected as optimum extraction time



Extraction at 5/10/15/20 minutes (n = 3)




Validation and Results

- Analyte distribution


Sintered (%)	CFP	Down. Tube	Front Imp.	Bridge	Back Imp.
Limonene	2.4	3.9	93.0	<LOD	0.7
Neomenthol	87.8	0.8	11.4	<LOD	<LOD
Menthol	95.1	0.5	4.4	<LOD	<LOD

Non-Sintered (%)	Front Imp. + CFP	Back Imp.
Limonene	93.8	6.2
Neomenthol	98.9	1.1
Menthol	99.6	0.4

- Non-sintered - up to ~6% carry-over to back impingers.



Non-Sintered



Sintered

Validation and Results

- Assessing Transfer (n = 5)

Sintered Impingers	Mean Transfer (%)	SD Transfer	%RSD	U
Limonene	50	4	9	9
Neomenthol	87	5	6	10
Menthol	93	5	5	10
Non-Sintered Impingers	Mean Transfer (%)	SD Transfer	%RSD	U
Limonene	54	2	3	3
Neomenthol	92	2	2	4
Menthol	97	2	2	4
Quantitative Results	Mean Transfer (%)	SD Transfer	%RSD	U
Limonene	63	6	9	11
Neomenthol	82	1	2	2
Menthol	87	1	1	2

Validation and Results

- Assessing Transfer (n = 5)

Sintered Impingers	Mean Transfer (%)	SD Transfer	%RSD	U
Limonene	50	4	9	9
Neomenthol	87	5	6	10
Menthol	93	5	5	10
Non-Sintered Impingers	Mean Transfer (%)	SD Transfer	%RSD	U
Limonene	54	2	3	3
Neomenthol	92	2	2	4
Menthol	97	2	2	4
Quantitative Results	Mean Transfer (%)	SD Transfer	%RSD	U
Limonene	63	6	9	11
Neomenthol	82	1	2	2
Menthol	87	1	1	2

Validation and Results

- Assessing Transfer (n = 5)

Sintered Impingers	Mean Transfer (%)	SD Transfer	%RSD	U
Limonene	50	4	9	9
Neomenthol	87	5	6	10
Menthol	93	5	5	10
Non-Sintered Impingers	Mean Transfer (%)	SD Transfer	%RSD	U
Limonene	54	2	3	3
Neomenthol	92	2	2	4
Menthol	97	2	2	4
Quantitative Results	Mean Transfer (%)	SD Transfer	%RSD	U
Limonene	63	6	9	11
Neomenthol	82	1	2	2
Menthol	87	1	1	2

Validation and Results

Sintered vs Non-Sintered vs Quantitative

- All methods show %RSD of $\leq 10\%$ (n = 5).
- Comparable results to quantitative method.
- Non-sintered impingers showed up to 6% break-through to back impinger. However:
 - Non-sintered - Mean % transfers were higher.
 - Non-sintered - Best overall %RSDs.
 - Non-sintered - Much easier to clean!
- More variability associated with sintered impingers – this section is hand made.

Validation and Results

- Injection precision


- Injection repeatability acceptable for both measured peak area and PAR.

Compound	Peak Area %RSD	PAR %RSD
Limonene	2.60	0.34
Neomenthol	2.08	1.30
Menthol	2.29	1.43
Anisole-D8 (Int. Std.)	2.48	N/A

- GC/MS sample carry-over in solvent blanks < 1% assessed by peak area.

Example Aerosol Chromatogram



- Target compounds - selected ion monitoring (SIM)
- Stabilwax-DA: 30m x 0.25mm x 0.5 μ m GC column




The Limit Test

- Are target compounds above or below defined thresholds (Limit of Concern) taking into account the uncertainties?

- Samples prepared so far:


- Aerosol impinger extracts (n = 5): 
- E-Liquid extracts (n = 5): 

- Limit Test sample preparation:



The Limit Test Red, Amber, Green (RAG)

- Possible outcomes:



Summary & Conclusions

- Presented a screening approach to flavour transfer:
 - CFP & liquid impingers
 - Use of peak area ratios for relative rather than absolute measurement of transfer

- Discussed concept of Limit Test:
 - RAG approach
 - Fail/succeed, fast/cheap

- The end of the calibration curve?
 - NO!

- Approach may also be applicable to tobacco heating products (THPs) and combustibles.

Acknowledgements and Reference

- Colleagues at BAT Analytical Development Centre

Stuart Martin



Malcolm Saxton



Louise Bishop



- Reference:

- [1] Costigan, S. & Meredith, C: An approach to ingredient screening and toxicological risk assessment of flavours in e-liquids, Regulatory Toxicology and Pharmacology, Volume 72, Issue 2, July 2015, Pages 361-369.

Thank you for listening, any questions?



Research &
Development

About

Science

News


Library

Careers



Discover BAT Science

An overview of our research activities



www.bat-science.com



@BAT_Sci

stuart_martin@bat.com

©British American Tobacco (Investments) Limited 2017. All rights reserved. No part of these materials may be reproduced in any form or by any means without the prior written consent of British American Tobacco (Investments) Limited and no responsibility or liability is accepted for any third party reliance on any data contained herein.

The data and information used in these materials has been compiled from a number of sources.