



Method Comparisons for Particle Size Distributions and Nicotine Dissolution Profiles in Smokeless Tobacco Products

Candice K. CUNNINGHAM*; Summer N. Hanna*; Christopher L. Keller#; Thomas R. Meadows#; Anzhi Kaszycki*

* - RAI Services Company (Winston Salem, NC, USA)

- RJ Reynolds Tobacco Company (Winston Salem, NC, USA)

72nd Tobacco Science Research Conference (Memphis, TN)

September 18, 2018

Particle Size & Dissolution in FDA Regulation

- **Pharmaceutical Industry**

- Particle size assessment and dissolution testing are recommended by FDA/CDER
- Generally recognized relationship between surface area and exposure to/release of “active” ingredient



- **Tobacco Industry**


- Based on pharma, FDA/CTP has hypothesized similar relationships exist for smokeless tobacco products (STPs)
- Currently, there is no empirical evidence in the scientific literature to support this hypothesis
- The nature of tobacco products, manner of use for STPs, and no claim of an “active” ingredient are strikingly dissimilar from pharmaceuticals



National/International Methods for Particle Size Assessment



Common Particle Size Assessment Methods



Laser Diffraction/Light Scattering

- Matrix challenges
- Distributions assessed typically assuming spherical particle shape




Sedimentation

- Matrix challenges
- Assumes spherical particle shape



Microscopy


- Representative sample challenges
- Distributions assessed typically assuming spherical or regular particle shape



Sieving


- Orientation through the screen
- Blinding

STP Particle Size Sieve Analysis Method



*U.S. Standard sieve sizes, ASTM Designation E11-17. "Standard Specifications for Woven Wire Test Sieve Cloth and Test Sieves." (2017)

Moist Snuff Product Particle Size Distributions



- Fine cut, long cut designation typically a tactile description
- Commercial Fine Cut STP distribution similar to CRP2
- Commercial Long Cut STP weight percent of sample shifted to the >1 mm particle size compared to Commercial Fine Cut STP

Comparing Particle Size Distributions

- K-S (Kolmogorov-Smirnov) two-sample test
- p-Value <0.05 indicates distributions are significantly different


Distributions Compared	p-Value
Commercial Fine Cut STP to Commercial Long Cut STP	<0.0001

National/International Methods for Dissolution




U.S. Food and Drug Administration Dissolution Methods: https://www.accessdata.fda.gov/scripts/cder/dissolution/dsp_getallData.cfm


Dissolution Apparatus for Determining Release Rates of Constituents



STP Nicotine Dissolution Method



Moist Snuff Product Dissolution



Comparing Dissolution Profiles

- Similarity factor (f_2): > 50 indicates “sameness”

$$f_2 = 50 \times \log \left[\frac{100}{\sqrt{1 + \frac{\sum_{t=1}^n (R_t - T_t)^2}{n}}} \right]$$

R_t = dissolution for Product 1 at time t

T_t = dissolution for Product 2 at time t



n = number of time points for which dissolution was assessed

Distributions Compared	f_2 Value
Commercial Fine Cut STP to Commercial Long Cut STP	65

Moore, JW, and Flanner, HH. “Mathematical Comparisons of Dissolution Profiles.” *Pharm. Tech.* June **1996**: 64 – 74.

Ma, M-C, Wang, BC, Liu, J-P, and Tsong, Y. “Assessment of Similarity Between Dissolution Profiles.” *J. Biopharma. Stat.* **2000**. 10 (2): 229 – 249.

Particle Size & Dissolution for Commercial STPs



Conclusions

- ❖ Fundamental approaches for comparing particle size distributions, and for dissolution profiles can be applied to evaluate smokeless tobacco products

- ❖ Method development requires specific considerations for smokeless tobacco products
 - ❖ Particle size
 - ❖ matrix effects
 - ❖ irregular particle shape
 - ❖ non-uniform and no “active” ingredient claim for assessment
 - ❖ Dissolution
 - ❖ samples do not disintegrate
 - ❖ numerous “ingredients” (i.e., tobacco constituents), but no “active” ingredient claim

- ❖ In these data presented, no differences were observed in nicotine release rate for products with statistically significantly different particle size distributions

Acknowledgments

Steven Flynn

Paul Ayres

Jeremi Johnson

Mike Borgerding

Mike Ogden

Walt Morgan

Summer Hanna

Hannah Carmichael

Cassandra West

ELL-PSS (Winston Salem, NC)

Catalent Pharma Solutions (Morrisville, NC)

90 232.03 Th Thorium	47 107.86 Ag Silver	7 14.00 N Nitrogen	19 39.09 K Potassium	16 32.06 S Sulfur
-----------------------------------	----------------------------------	---------------------------------	-----------------------------------	--------------------------------